
DENTASONIC

Products IFU

Instructions
for use
and safety
recommendation

ⓐ ⓑ ⓒ ⓓ ⓔ ⓕ ⓖ ⓗ ⓘ

Description

ErgoStrip are handheld strips that can be used to remove set cement, break contacts, Interproximal Reduction (IPR) and for polishing and finishing. ErgoStrip is available in straight and Resin diamond coated configurations (to better follow the contour of teeth). Choices include single-sided for small reduction, double-sided for larger reduction, and serrated for breaking contact and removing cement that has set in the interproximal space. They come in superfine, fine, medium and coarse grits.

- Diamond grit sizes: 15µm, 25µm, 40µm, 60µm
- Strip height: 3.6mm
- Length of the diamond-coated part: 20mm
- Sterilizable and reusable

Indications for Use

- ErgoStrip is indicated for the removal of set cement, to break contacts, for Interproximal Reduction (IPR) procedures, and for polishing and finishing.
- ErgoStrip (single-sided, double-sided and polyester) is indicated for both anterior and posterior teeth.
- Pass the strip back and forth gently and buccolingually only.

Contraindications

ErgoStrip is contraindicated if the Instructions for Use (IFU) is not strictly followed or with patients with poor enamel structure (hypoplasia), high caries index, poor oral hygiene or severe tooth crowding.

Warnings

- Avoid bending and levering: use of excessive force (strong levering) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Avoid temperatures above 180°C which may affect the durability of the instrument.
- Avoid applying forces greater than 3N as this could cause damage to the instrument.
- To prevent patient from ingesting or inhaling the ErgoStrip, prior to using ErgoStrip, thread a piece of dental floss through the arc of the ErgoStrip. Secure by tying a slipknot. Wrap the floss around one of your other fingers to secure the ErgoStrip.
- The consumption of the instrument lifetime until its worn out, replacement and disposal is an exclusive decision of the physician according to his operative needs.
- To ensure traceability during the instrument lifetime we recommend to preserve the packaging.

Principles of safety

- Contact pressure: bend the strip 1-1.5 mm to achieve the ideal 1-2.5N operating force.
- Do not stop the motion of the file while the same is in the interdental space
- In case of overcrowding open the interdental space using a wedge in the cervical area
- Do not bend, warp or leverage the instrument
- Not observing these precautions leads to a high risk in terms of safety and can lead to the possible breakage or damage of the instrument

Precautions

- Sterilize ErgoStrip prior to patient use. See details below.
- Read the Instructions for Use (IFU) completely before using. Use of ErgoStrip inconsistent with IFU is at the sole discretion and sole responsibility of the user.
- Do not use any ErgoStrip that is damaged, bent, broken, or show signs of wear.
- For maintaining anatomical contour in the posterior and to refine the natural anatomy of the anterior, finish with the curved version.
- Use a gap gauge frequently during the procedure when gap width accuracy is critical.

Adverse Reactions

- None known. However, ErgoStrip is made from stainless steel and contains nickel. Nickel may cause allergic reactions in some patients.

No Reactions for allergic patients

- Resin Diamond ErgoStrip is made from resin and **NOT** contains nickel. Resin will **NOT** cause allergic reactions to patients.

Instructions for Cleaning, Sterilization and Disposal

Warnings: DO NOT use cold sterile or dry heat sterilization processes.

Limitations on reprocessing: ErgoStrip can be sterilized up to 10 times. Resin ErgoStrip can be sterilized up to 3-5 times.

INSTRUCTIONS

Point of use: Sterilize before first use and between patients.

Preparation for decontamination: It is best to wipe immediately after use to avoid contaminants drying, thereby making the cleaning process more time-consuming and difficult. Immerse or soak ErgoStrip in a hospital grade enzymatic detergent such as enzymatic solution for 15 minutes.

Cleaning: Manual Clean thoroughly - Scrub all external surfaces with a soft bristle brush until all visible soil has been removed. Use a small diameter brush cleaner to clean a cannulation holes. Inspect for visible soil on exposed surfaces. Ultrasonically clean using a solution that is recommended for dental instruments such as an enzymatic solution for 5 – 10 minutes. Rinse ErgoStrip under copious amount of warm water for 2 -3 minutes and visually inspect to confirm removal of debris. Do not use highly aggressive chemical products (e.g. hydrochloric acid, hydrogen peroxide), as they could corrode the instruments.

These substances can also impair the technical properties of the plastics, i.e. change their hardness and durability.

Cleaning: Automated Recommended

Disinfection: Recommended

Drying: Allow to air dry in clean area. Blow lumens with clean air using filtered air source or syringe.

Maintenance, inspection and testing: Visually inspect for fractured or worn abrasive strip that could prevent proper function. Discard in Sharps container.

Packaging: Bag in an FDA cleared sterilization pouch, Sterilization pouches are recommended.

Sterilization: EU and US Autoclave Parameters

Steam autoclave- Gravity (Wrapped)

Minimum Exposure Temperature: 121°C (250°F)

Minimum Exposure Time: 30 minutes

Minimum Dry: 30 minutes

EU Autoclave Parameters

Steam autoclave- Gravity (Unwrapped)

Minimum Exposure Temperature: 134°C

Minimum Exposure Time: 3 minutes

Minimum Dry: None

US Autoclave Parameters

Steam autoclave- Pre-Vacuum (Wrapped)

Minimum Exposure Temperature: 135°C (275°F)

Minimum Exposure Time: 3 minutes

Minimum Dry: 30 minutes

Storage: ErgoStrip should be stored in the sterilization pouch until ready for use with a maximum storage period of 6 months before product needs to be re-sterilized. Unwrapped ErgoStrip should be used immediately after sterilization.

Instructions for Use
Cement Removal and Break Contacts

Use ErgoStrip Saw (Blue) serrated strips to remove set cement and break contacts.

- Hold the ErgoStrip between your thumb and index finger.
- Insert the serrated blade between the teeth and with a gentle "sawing motion", to remove the cement or to break contacts.

Finishing and Polishing

Use ErgoStrip Single-Sided Abrasive Strips for both anterior and posterior teeth. When following the contour of teeth is critical, use the Resin Abrasive Strips as they are better suited to this purpose.

- Select the desired grit metal strip (White – SuperFine, Yellow – Fine).
- Select the desired grit Resin Strips (Blue – UltraFine, Yellow – Superfine, Green – Fine, White – Medium or Gray – Coarse).
- Hold the ErgoStrip or the Resin ErgoStrip between your thumb and index finger.
- In an angled position, begin a gentle "sawing motion" until light resistance is felt.
- If needed, continue to other grits until the proper finish/polish is achieved.

Interproximal Reduction (IPR)

Use Single-sided and Double-Sided ErgoStrip in the sequence shown below for IPR procedures. When following the contour of teeth is critical, use the Resin ErgoStrip Strips as they are better suited to this purpose. For creating IPR spaces from 0.08mm to 0.2mm, use the Single-Sided ErgoStrip Abrasives; for IPR spaces greater than 0.2mm, use the Double-Sided ErgoStrip Abrasives.

- Regardless of the strip type chosen (single, double-sided, straight metal or Resin) as described above, begin with the superfine grit (White).
- Hold between the thumb and index finger. In an angled position, begin a gentle "sawing motion" until very light resistance is felt.
- It is best to check the gap with a verified gap measuring gauge before moving to the next color (grit) in the sequence. Once the desired gap is achieved there is no need to continue to the next larger grit size.
- If additional space is needed, switch to the fine grit (Yellow) and use in the same manner as described above. Check gap width.
- If additional space is needed, switch to medium grit (Red). Check gap width.
- If additional space is needed, switch to the coarse grit (Blue). Check gap width. Continue until the desired gap is achieved.
- If needed, continue to lesser grits until the proper finish/polish is achieved.

IPR Sequence Guide For Orthodontic Procedures

In order to create precise interproximal spaces for Orthodontic procedures, use ErgoStrip color coded strips according to the sequence outlined below. The reduction process (IPR) should begin with the white single-sided grit (0.08 mm), and continue in sequence until the desired space is achieved. Check the reduction with an approved gap measurement gauge before moving to the next color (grit) in the sequence. If needed, continue to lesser grits until the proper finish/polish is achieved.

Note: Use all versions of ErgoStrip (single-sided – ESM1, double-sided – ESM, and Resin – ES) in any area of the mouth. For maintaining anatomical contour in the posterior and to refine the natural anatomy of the anterior, finish by using the Resin Diamond coated ErgoStrip version.

It is necessary to report any serious accident occurring in relation to the medical device supplied by us to the manufacturer and to the competent authority of the Member State in which it is based.

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DentaSonic ErgoStrip, Ref. ESM-15
 DentaSonic ErgoStrip, Ref. ESM-25
 DentaSonic ErgoStrip, Ref. ESM-40
 DentaSonic ErgoStrip, Ref. ESM-60

DentaSonic ErgoStrip, Ref. ESM1-15
 DentaSonic ErgoStrip, Ref. ESM1-25
 DentaSonic ErgoStrip, Ref. ESM1-40
 DentaSonic ErgoStrip, Ref. ESM1-60

DentaSonic ErgoBee, Ref. ESB-15
 DentaSonic ErgoBee, Ref. ESB-25
 DentaSonic ErgoBee, Ref. ESB-40
 DentaSonic ErgoBee, Ref. ESB-60

DentaSonic ErgoSaw, Ref. ESS-15
 DentaSonic ErgoSaw, Ref. ESS-25
 DentaSonic ErgoSaw, Ref. ESS-40

DentaSonic ErgoSaw, Ref. ESS-1
 DentaSonic ErgoPlaque, Ref. ESP-1

DentaSonic Resin ErgoStrip, Ref. ES-9
 DentaSonic Resin ErgoStrip, Ref. ES-15
 DentaSonic Resin ErgoStrip, Ref. ES-30
 DentaSonic Resin ErgoStrip, Ref. ES-45
 DentaSonic Resin ErgoStrip, Ref. ES-60

DentaSonic Resin CarbideStrip, Ref. RECS-15

	Simboli etichette DentaSonic	DS Produktetikette Symbole	DS Symbol des etiquettes	DS Product symbol
	Fabbricante	Produzent	Fabricant	Manufacturer
	Monouso	Einweg Produkt	Jetable	Disposable
	Non sterile	Nicht Steril	Pas sterile	Non sterile
	Attenzione	Achtung	Attention	Attention
	Informazioni	Informationen	Informations	Information
	Riciclaggio	Recycling	Recyclage	Recycling
	Smaltimento controllato	Kontrollierte Entsorgung	Elimination contrôlée	Controlled elimination
	Numero di Lotto	LOT Nummer	Numero de LOT	LOT number
	Rappresentante Europeo	EU Vertreter	Représentant Européen	EU Representative

Packaging of Kit
Kit-Inhalt
Emballages par Set
Confezioni da Set
Embalajes de kit

DentaSonic ErgoStrip, Set, Ref. DSM-55
 (2xDSM-15, DSM-25, DSM-40, DSM-60)

DentaSonic ErgoStrip, Set, Ref. DSM1-55
 (2xDSM1-15, DSM1-25, DSM1-40, DSM1-60)

DentaSonic ErgoSaw Set 5, Ref. DSS-55
 (2xDSS-15, 2xDSS-25, 1xDSS-40)

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Quaser S.r.l.