

WE
CARE FOR
ENAMEL



DENTA
SONIC

Packaging refills, set and Kit's for your orders



DiaStrip 2 sides and 1 side.



Name:	DiaStrip	DiaStrip	DiaStrip	DiaStrip	DiaStrip
PN: 2 sides:	DSM-01-12	DSM-02-12	DSM-03-12	DSM-04-12	RDS-05-12
PN: 1 side left:	DSM-01L-12	DSM-02L-12	DSM-03L-12		
PN: 1 side right:	DSM-01R-12	DSM-02R-12	DSM-03R-12		
Thickness mm:	0.1mm	0.2mm	0.3mm	0.4mm	0.5mm

Diamond Resin DiaStrip 2 sides and 1 side.



Name:	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip
PN: 2 sides:	RDS-02-12	RDS-025-12	RDS-03-12	RDS-04-12	RDS-05-12
Thickness mm:	0.2mm	0.25mm	0.3mm	0.4mm	0.5mm
PN: 1 side left:	RDS-008L-12	RDS-009L-12	RDS-010L-12	RDS-012L-12	
PN: 1 side right:	RDS-008R-12	RDS-009R-12	RDS-010R-12	RDS-012R-12	
Thickness mm:	0.08mm	0.09mm	0.10mm	0.12mm	

DiaSaw 2 sides and 1 side.



Name:	DiaSaw	DiaSaw	DiaSaw
PN: 2 sides:	DSS-01-12	DSS-02-12	DSS-1-12
Thickness mm:	0.01mm	0.02mm	0.05mm
PN: 1 side left:	DSS-005L-12		not coated
PN: 1 side right:	DSS-005R-12		
Thickness mm:	0.05mm		

DiaSaw no Diamond.



Name:	DiaSaw
PN: 2 sides:	DSS-1-12
Thickness mm:	0.05mm
PN: 1 side left:	not coated
PN: 1 side right:	
Thickness mm:	

DiaBee 2 sides.



Name:	DiaBee	DiaBee	DiaBee	DiaBee	DiaPlaque
PN: 2 sides:	DSB-15-12	DSB-25-12	DSB-40-12	DSB-60-12	DP-1-12
Grain µm:	15µm	25µm	40µm	60µm	not coated

Diaplaque no Diamond.



Name:	Diaplaque
PN: 2 sides:	DP-1-12
Grain µm:	not coated

To order blister of only 1 or 3 strips, change the end of the code, from -12 to -1 or -3

DiaStrip 2 sides and 1 side.



Name:	DiaStrip	DiaStrip	DiaStrip	DiaStrip	DiaStrip
PN: 2 sides:	DSM-01	DSM-02	DSM-03	DSM-04	RDS-05
PN: 1 side left:	DSM-01L	DSM-02L	DSM-03L		
PN: 1 side right:	DSM-01R	DSM-02R	DSM-03R		
Thickness mm:	0.1mm	0.2mm	0.3mm	0.4mm	0.5mm

Diamond Resin DiaStrip 2 sides and 1 side.



Name:	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip
PN: 2 sides:	RDS-02	RDS-025	RDS-03	RDS-04	RDS-05
Thickness mm:	0.2mm	0.25mm	0.3mm	0.4mm	0.5mm
PN: 1 side left:	RDS-008L	RDS-009L	RDS-010L	RDS-012L	
PN: 1 side right:	RDS-008R	RDS-009R	RDS-010R	RDS-012R	
Thickness mm:	0.08mm	0.09mm	0.10mm	0.12mm	

DiaSaw 2 sides and 1 side.



Name:	DiaSaw	DiaSaw	DiaSaw
PN: 2 sides:	DSS-01	DSS-02	DSS-1
Thickness mm:	0.01mm	0.02mm	0.05mm
PN: 1 side left:	DSS-0.05L		not coated
PN: 1 side right:	DSS-0.05R		
Thickness mm:	0.05mm		

DiaSaw no Diamond.



Name:	DiaSaw
PN: 2 sides:	DSS-1
Thickness mm:	0.05mm
PN: 1 side left:	not coated
PN: 1 side right:	
Thickness mm:	

DiaBee 2 sides.



Name:	DiaBee	DiaBee	DiaBee	DiaBee	Diaplaque
PN: 2 sides:	DSB-15	DSB-25	DSB-40	DSB-60	DP-1
Grain µm:	15µm	25µm	40µm	60µm	not coated

Diaplaque no Diamond.



Name:	Diaplaque
PN: 2 sides:	DP-1
Grain µm:	not coated

PoliStrip 2 sides and 1 side. Final polishing






Name:	PoliStrip	PoliStrip
PN: 2 sides:	RPS-9	RPS-15
PN: 1 side left:	RPS-9L	RPS-15L
PN: 1 side right:	RPS-9R	RPS-15R
Grain µm:	9µm	15µm
Content:	12X RPS-9	12X RPS-15

PoliStrip 2 sides and 1 side. Final polishing



Name:	PoliStrip	PoliStrip
PN: 2 sides:	RPS-S2	RPS1-S2
Grain µm:	9µm, 15µm	9µm, 15µm
Content:	6X RPS-9, 15	6X L or R RPS-9, 15

Polishing enamel surface result



Untreated enamel surface and treated surface with Resin PoliStrip 15µm SEM enlargement 500x

Untreated enamel surface and treated with Resin PoliStrip 15µm SEM enlargement 2000x


Treated enamel surface (40µm) and center treated with Resin PoliStrip 15µm SEM enlargement 2000x



Courtesy of the Institute of Material Science: Technische Universität Darmstadt, Germany

How insert and remove correctly the strip from the handpiece

Wrong strip handling



Never push inside the strip into the head like the two picture on the left, because you will not have enough force, and you risk to break the shank /blade of the strip.

Correct strip handling



To insert the strip proceed by push the shank into the head, by turning on the right and left side until the grip are ok.
To remove the strip proceed by pull out from the head by turning on the right and left side.
 Never plie the strip because you can broke the shank. The shank resist until 1kg force by plie it.

DiaStrip 2 sides and 1 side.



Name:	DiaStrip	DiaStrip	DiaStrip	DiaStrip
PN: 1/2 sides:	DSM-S12	DSM-S55	DSM1-S12	DSM1-S55
Content:	2X DSM-0.1, 0.2, 0.3, 0.4 2X RDS-0.5, 2X RPS-15	1X DSM-0.1, 0.2, 0.3, 0.4, 1X RDS-0.5	2X DSM-0.1L, 0.1R, 0.2L, 0.2R, 0.3L, 0.3R	1X L or R DSM-0.1, 0.2, 0.3 2X RPS-15

Resin DiaStrip 2 sides.



Name:	Resin DiaStrip
PN: 1/2 sides:	RDS-S12
Content:	2X RDS-0.2, 0.25, 0.3, 0.4, 0.5 1X RPS-9, 15

DiaSaw 2 sides and 1 side.



Name:	DiaSaw	DiaSaw
PN: 1/2 sides:	DSS-S12	DSS-S55
Content:	1X DSS-0.05L, 0.05R, 4X 0.1, 2X 0.2, 2X RPS-9, 15	1X DSS-0.05L, 0.05R, 0.1, 0.2, 1X RPS-15

DiaBee 2 sides.



Name:	DiaBee
PN:	DSB-S55
Content:	1X DSB-15, 25, 40, 60 1X RPS-15

Aligner DiaStrip 2 sides.



Name:	Aligner DiaStrip	Aligner DiaStrip
PN: 1/2 sides:	DSA-S12A	DSA-S5A
Content:	2X DSS-0.1, 2X DSM-0.2, 0.3 1X DSM-0.4, 0.5 1X RDS-0.4, 0.5 2X RPS-15	1X DSS-0.1, 1X DSM-0.2, 0.3, 0.4 1X RDS-0.5

ErgoStrip Packaging & refills - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side



Name:	ErgoStrip	ErgoStrip	ErgoStrip	ErgoStrip
PN: 2 sides:	ESM-15-12	ESM-25-12	ESM-40-12	ESM-60-12
Content:	12X ESM-15	12X ESM-25	12X ESM-40	12X ESM-60
PN: 1 side:	ESM1-15-12	ESM1-25-12	ESM1-40-12	ESM1-60-12
Content:	12X ESM1-15	12X ESM1-25	12X ESM1-40	12X ESM1-60
Grain µm:	15µm	25µm	40µm	60µm

ErgoSaw 2 sides.



Name:	ErgoSaw	ErgoSaw	ErgoSaw	ErgoSaw
PN: 2 sides:	ESS-15-12	ESS-25-12	ESS-40-12	ESS-1
Content:	12X ESS-15	12X ESS-25	12X ESS-40	12X ESS-1
Grain µm:	15µm	25µm	40µm	not coated

ErgoBee 2 sides.



Name:	ErgoBee	ErgoBee	ErgoBee	ErgoBee
PN: 2 sides:	ESB-15-12	ESB-25-12	ESB-40-12	ESB-60-12
Content:	12X ESB-15	12X ESB-25	12X ESB-40	12X ESB-60
Grain µm:	15µm	25µm	40µm	60µm

until end of stock

ErgoPlaque no Diamond.



Name:	ErgoPlaque
PN: 1/2 sides:	EP-1
Content:	12X EP-1
Grain µm:	not coated



Ergostrip Refill & Set confection.

ErgoStrip Packaging & Set's - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side.



Name:	ErgoStrip	ErgoStrip	ErgoStrip
PN: 1/2 sides:	ESM-S12	ESM1-S12	ESM-S21
Content:	3X ESM-15, 25, 40, 60	3X ESM1-15, 25, 40, 60	1X ESM1-15, 25, 40, 60 2X ESM-15, 25, 40, 60

ErgoSaw 2 sides.



Name:	ErgoSaw
PN: 1/2 sides:	ESS-S12
Content:	4X ESS-15, 25, 40

ErgoBee 2 sides.



Name:	ErgoBee
PN: 1/2 sides:	ESB-S12
Content:	3X ESB-15, 25, 40, 60 <i>until end of stock</i>

Packaging & Refills - Restoration, Veneers, Crown & Bridges

Polishing ErgoStrip 2 sides, and 1 side.



Name:	Polishing ErgoStrip
PN: 1/2 sides:	REPS-15
Content:	12X REPS-15
Grain µm:	15µm

until end of stock

How use ErgoStrip



Finger Grip design for better tactile control & perfect Grip.

ErgoStrip Packaging & Refills - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side



Name:	ErgoStrip	ErgoStrip	ErgoStrip	ErgoStrip
PN: 2 sides:	ESM-15-4	ESM-25-4	ESM-40-4	ESM-60-4
Content:	4X ESM-15	4X ESM-25	4X ESM-40	4X ESM-60
PN: 1 side:	ESM1-15-4	ESM1-25-4	ESM1-40-4	ESM1-60-4
Content:	4X ESM1-15	4X ESM1-25	4X ESM1-40	4X ESM1-60
Grain µm:	15µm	25µm	40µm	60µm

ErgoSaw 2 sides.



Name:	ErgoSaw	ErgoSaw	ErgoSaw	ErgoSaw
PN: 2 sides:	ESS-15-4	ESS-25-4	ESS-40-4	ESS-1-4
Content:	4X ESS-15	4X ESS-25	4X ESS-40	4X ESS-1
Grain µm:	15µm	25µm	40µm	not coated

ErgoBee 2 sides.



Name:	ErgoBee	ErgoBee	ErgoBee	ErgoBee
PN: 2 sides:	ESB-15-4	ESB-25-4	ESB-40-4	ESB-60-4
Content:	4X ESB-15	12X ESB-25	4X ESB-40	4X ESB-60
Grain µm:	15µm	25µm	40µm	60µm

until end of stock

ErgoStrip Packaging & Set's - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side.



Name:	ErgoStrip	ErgoStrip
PN: 1/2 sides:	ESM-S4	ESM1-S4
Content:	1X ESM-15, 25, 40, 60	1X ESM1-15, 25, 40, 60

ErgoSaw 2 sides.



Name:	ErgoSaw
PN: 1/2 sides:	ESS-S4
Content:	1X ESS-15, 25, 40, 1

ErgoBee 2 sides.



Name:	ErgoBee
PN: 1/2 sides:	ESB-S4
Content:	1X ESB-15, 25, 40, 60

until end of stock

Basic Set's

Basic Set Ortho

Art. Ref.:	DSG4-204
Stroke:	4:1
Pieces:	14
Sterilizable:	Yes
Application:	Blue, Pink

with a 1 additional reciprocating head
Warranty option: 1 or 2 years



1 additional reciprocating head

Silent handpiece, no vibration. Internal water cooling and interchangeable head

Basic Set Aligner

Art. Ref.:	DSG4A-204
Stroke:	4:1
Pieces:	13
Sterilizable:	Yes
Application:	Blue, Pink

Warranty option: 1 or 2 years



Basic Set Orthofamily

Art. Ref.:	DSG4-204-OF
Stroke:	4:1
Pieces:	13
Sterilizable:	Yes
Application:	Blue, Pink

Warranty option: 1 year



Tested by Tor Vergata University, Rome.

Recognised Invisalign center for Invisalign certification courses.

Basic Set Ortho NWC

Art. Ref.:	DSG4-200
Stroke:	4:1
Pieces:	12
Sterilizable:	Yes
Application:	Blue, Pink

Warranty option: 1 year



No Water-cooled Handpiece

Professional Kit 1

Art. Ref.: DSPK1-2Y
 Stroke: 4:1
 Pieces: 35
 Sterilizable: Yes
 with 1 additional rotating push button head
 Warranty option: 2 years



Professional Kit 2



Art. Ref.:	DSPK2-2Y	DSPK2-1Y	DSPK2L-2Y
Stroke:	4:1	4:1	4:1
Pieces:	35	35	35
Sterilizable:	Yes	Yes	Yes
	with 1 additional reciprocating head, and voucher	no voucher	with voucher
Warranty option:	2 years	1 year	2 years



Professional Kit 3

Art. Ref.: DSPK3-2Y
 Stroke: 4:1
 Pieces: 36
 Sterilizable: Yes
 with 1 additional reciprocating and 1 rotating head
 Warranty option: 2 years



Take
 the guesswork
 out of
IPC interproximal
 correction



DENTASONIC

For more information about all this products, please contact us online:
5% off for all online order with code "DS5" on www.denta-sonic.com

DS DENTASONIC SA
Via Piodella 12
CH-6933 Muzzano
Switzerland

Tel. +41 91 600 19 60
info@dentasonic.ch
www.dentasonic.ch

Official Distributor:

European Representative: QuaSer S.r.l., Via Pian due Torri, 19, 00146 Roma

Some of the products and designations mentioned in the text are patented or copyrighted. The absence of a special reference or the ® sign should not be interpreted as the absence of copyright protection. This publication is copyrighted. All rights, including those of translation, reprinting, and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm, or other methods) without the editor's prior consent in writing. Colours and products subject to alterations. Printing errors excepted.